

Vikan - ultra microfibre technology



Vikan Original cloth,
cleaning cloth

Article	Article no.	Size	Recommended frame or holder
Green	691012	310x320 mm	
Blue	691013	310x320 mm	
Red	691014	310x320 mm	
Yellow	691016	310x320 mm	
Traceability	First 4 numbers describes the construction and the last numbers indicates the model of the product. Batch number indicates when the product is manufactured, referring to this Vikan can trace the production and take action if necessary.		
Packed	5 pcs. in a plastic bag - 50 pcs. in a carton.		

CLEANING

Where to use:	Cleaning of all types of surfaces with high demands on cleaning efficiency and hygienic. Perfect for regular maintenance cleaning. Avoid any contact with solvents, strong acid or alkaline products.		
How to use:	Preferably used in Vikan Vikan® system. The cloth can be used either damp or dry - wipe with the cloth folded - fold in another direction until every side has been used as face side - thereafter released in the laundry bag. Training is organized by Vikan. Contact our sales organization for further information or set up a meeting.		
Optimum moisture:	75% to achieve best cleaning efficiency and scope		
Preparation:	Dry or damp. Be aware of the recommended optimum moisture.		
Article	Dry weight ± 5 g	ml water	Damp weight maximum
Original cloth	39	29	68
Drying speed	The controlled adding of damp results in an uniform and quick drying surface. The cleaning result is very obvious and immediate. Drops or stripes can be avoided.		

CARE

Care label			
Care label describes the maximum temperature that should be used. Due to cost, time and environmental impact a reduced temperature should always be evaluated depending on the application.			
Recommended washing	70°C	Use detergents developed for microfiber products with a pH less than 10. Softner or biocides should not be used as they can cause fiber damage and specific in combination with high temperature. The electrostatic characteristic that is essential for the dust collecting will be reduced or lost with this type of treatment.	
Washing in specific hygiene use	72-95°C, depending on time		
Drying instruction	Max 55°C	The recommendation is made to avoid the textile to be damaged if humidity is not controlled.	
Autoclaveable	steam at 121° C 20 min	Possible to autoclave in steam if required. Other methods must be tested and approved by the originator.	
Disposal after use	Material recovery	The product is pure synthetic with high strength and can be used in many re-circulated materials depending on their specification.	
	Energy recovery	The fibers have a high energy value and can be totally burned without anything else than CO ₂ and water.	



Vikan - ultra microfibre technology

COMPOSITION

Cleaning surface	80 % PES / 20 % PA. Splitted micro fiber: <0,3dtex, ACTEX®.	Uniform and dense loops of microfiber provides large contact area between the cloth and cleaning surface. The more contact the better efficiency in the cleaning. Our loops creates a large volume and helps to collect dust and dirt effecient. Fabric is provided with Vikan wicking performance which permanently protect the fibre and increase the transport inside the mop and enhance the cleaning and laundry efficiency.
------------------	---	---

CHARACTERISTICS

Shrinkage	Less than 5% at recommended washing conditions.	
Absorbency	Fast and immediate wicking properties. Can easily keep 5x its own weight in moisture.	
Staining	No staining in dry or damp condition if recommendation is followed. Some staining can be noticed on other textiles when the higher temperature range is used but separate washing of mops and cloths that are of different construction or used for different applications are recommended.	
Eco Label	Nordic Eco label for Micro fibre cloths and mops.	

GUARANTEE AND PERFORMANCE

500 times wash and use guarantee in general or hygiene use. Depending on cleaning areas (type of dirt and frequency of cleaning) the guarantee must be limited - but selecting the right mop and care, the life cycle can be optimized to even higher level. Ask our experts for best choice.

Lenience	Our micro fibre construction is tested against several glossy surfaces to secure the leniency of the mops and cloths. But before use on specific sensitive surfaces check that no particles as sand or old dirt are kept in the textile. Use a CD or similar as a test surface if you are not convinced.	
<i>Weber&Leucht S301/1.3</i>	Test is carried out at a external laboratory in order to meet the Nordic Ecolabel requirements on cleaning efficiency. In average the cloth recieved 25% better cleaning performance compared with standard methods applied with common chemicals. The analyzed product is ecological as it does not require cleaning agents - while still offering a distinctly enhanced cleaning efficiency.	
Bacteria reduction	Several independant test have been carried out with the common result that Vikan Original cloth is the absolut best in class to reduce bacteria and dirt on different surfaces. The ultra fine micro fiber has the exactly right size for attacking bacteria and fat, dirt and other nutrients where they can multiply.	

PRODUCT FEATURES AND BENEFITS

Ergonomics	Low weight - folded for better handle and utilization of the cloth surface	
Efficiency	High hygiene level - efficient reduction of bacteria - no spreading out of dirt - fast drying without residual chemicals	
Effective and simplicity	Easy to use - no changing of water or detergent - color coding system	
Environmental impact	No chemicals - high durability - effective utilized resources in manufacturing and use	



Developed and manufactured by the Vikan group for professionals.

