

# AURELIA® Delight PF®

Examination Gloves



-	(XS)	ø 100	10 x 100
<b>38226</b>	(S)	ø 100	10 x 100
<b>38227</b>	(M)	ø 100	10 x 100
<b>38228</b>	(L)	ø 100	10 x 100
<b>38229</b>	(XL)	ø 100	10 x 100



Vinyl



Powder Free



Clear



5.0 +/- 0.2 g



Smooth




3 Years

The gloves are made from a very soft, thin vinyl which allows for good, tactile sensitivity for all applications.

Beaded cuff, Ambidextrous, Glueless, Free of Thiurams  
Straight fingers, Non-Sterile

## Sizes & Dimensions



<b>XS</b>	75 ± 10 mm	min. 240 mm
<b>S</b>	80 ± 10 mm	min. 240 mm
<b>M</b>	95 ± 10 mm	min. 240 mm
<b>L</b>	110 ± 10 mm	min. 240 mm
<b>XL</b>	≥ 110 mm	min. 240 mm

## Single wall thickness



- ① min. 0,06 mm
- ② min. 0,06 mm
- ③ min. 0,05 mm

## Physical properties

	Unaged	Aged
Tensile Strength	min. 11,0 MPa	min. 11,0 MPa
Elongation	min. 350 %	min. 300 %
Strength at Break	min. 3,6 N	min. 3,6 N

## Quality Signs & Certificates

89/686/EEC, 93/42/EWG

ISO 9001:2008, ISO 13485, CMCAS 13485 (SGS)

EN 455 1-4

EN 420:2003 + A 1:2009

EN 1186: Food Contact Approved (SGS)

