

PRODUCT SPECIFICATION

Product: SALVEQUICK FABRIC (TEXTILE)
ADHESIVE BANDAGE (PLASTER)

Reference code: 6444CAP



1. General description

Adhesive bandages of two sizes, medium and senior, packed in a cardboard refill. The sterile adhesive bandages are covered by a paper cap. The refill fits into the Cederroth (Salvequick) dispenser for adhesive bandages.

The adhesive bandage consists of a fabric backing, evenly spread with an adhesive coating. An absorbent pad is fixed centrally to the dressing leaving a margin of adhesive surface around the pad. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When taking the adhesive bandage out from the package, the release-coated plastic foil will remain behind. The cap is still on the dressing. One half of the dressing can now be applied on the skin using one hand only. Then the cap and the rest of the coated paper can be removed and the adhesive bandage is ready for use.

2. Intended use

The sterile adhesive bandages are intended to be used as a mechanical barrier and for absorption of exudates from minor wounds.

3. Technical description

Product dimensions

Medium	Adhesive bandage:	72x19 mm
	Wound pad:	24x12mm
Senior	Adhesive bandage:	72x25 mm
	Wound pad:	24x17 mm

Material

Backing	Elastic, skin coloured, textile fabric
Adhesive	Hypoallergenic hotmelt adhesive
Woundpad	Dry laid, needle punched non woven covered with a wound contact layer
Facing paper	Single sided siliconised release coated paper/foil
Wrapping paper	Plain coldseal paper + paper cap

4. Sterilization

Sterilization method: Gamma irradiation

5. Composition

Refill

<u>Ref #</u>	<u>Product/Material</u>	<u>Pcs</u>
45165020	Medium Fabric adhesive bandage with cap	24
45165030	Senior Fabric adhesive bandage with cap	16
Dimensions (l x w x h) mm:		30 x 87 x 100

Consumer package

Number of refills in consumer package:	6
Dimensions (l x w x h) mm:	100 x 150 x 90

Transport package

Number of consumer packages in transport package:	10
Dimensions (l x w x h) mm:	520 x 160 x 190
Weight (kg):	2,8

6. Labeling

Product (refill)

Cederroth address, telephone- and faxnumber
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Symbols: CE +NB#, STERILE, R, DO NOT REUSE, LOT, USE BY

Consumer package

Cederroth address, telephone- and fax-number, web- and mail-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Quantity
 Symbols: CE +NB#, LOT, USE BY

Transport package

Cederroth address, telephone- and faxnumber, web-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Quantity

7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

8. Accessories

The following accessories can be used with the product;

Ref #: 1907 Salvequick Plaster Dispenser

Authorized by: Teija Ålander
Quality Assurance Manager

Date: 2008-02-01